


### 3 The Preparation of Buffers at Desired pH




**Objectives:** To become familiar with operating a pH meter, and to learn how to use the Henderson-Hasselbalch equation to make buffer solutions at a desired pH value.

**Introduction:** A buffer system is a mixture of a weak acid or a weak base and its salt (conjugate base or conjugate acid, respectively) that permits solutions to resist large changes in pH upon the addition of small amounts of hydrogen ions ( $H^+$ ) or hydroxide ions ( $OH^-$ ). What does this mean? It means a buffer helps maintain a *near constant pH* upon the addition of small amounts of  $H^+$  or  $OH^-$  to a solution.


In the definition of a buffer system, it should be noted that reference is made to weak acids and bases. What is the distinction between a *weak* and a *strong* acid or base? **Strong** acids (e.g., HCl) and bases (e.g., NaOH) almost completely ionize in water:




Many acids and bases, however, do not undergo complete ionization in water. These compounds, then, are called **weak** acids and bases. Many organic acids are weak acids. Weak acids do not completely dissociate in water. The dissociation of an organic acid is described by the following reaction:



Note that the unprotonated product of the acid dissociation reaction is referred to as a conjugate base. Recall the Bronsted-Lowry definition of acids and bases. An acid is defined as a proton ( $H^+$ ) donor, whereas the base is the  $H^+$  acceptor. In a conjugate acid/base pair such as that indicated above, the weak acid always has *one more*  $H^+$  than its conjugate base. Likewise, a weak base will have *one less*  $H^+$  than its conjugate acid.





One can analyze the strength of a weak acid. This means that the amount of hydrogen ion released can be determined. To do this, one can use the following expression:

$$K_a = \frac{[\text{H}^+][\text{A}^-]}{[\text{HA}]}$$

where  $K_a$  is defined as the acid dissociation constant.

The larger the value of  $K_a$ , the stronger the acid is. Because  $K_a$  values vary over a wide range, they are usually expressed using a logarithmic scale:

$$\text{pKa} = -\log K_a$$

The hydrogen ion is one of the most important ions in biological systems. The concentration of this ion affects most cellular processes. For example, the structure and function of most biological macromolecules and the rates of most biochemical reactions are strongly affected by  $[\text{H}^+]$ . The pH scale has been devised as a convenient method of expressing hydrogen ion concentration. pH has been defined as the negative logarithm of the hydrogen ion concentration:

$$\text{pH} = -\log [\text{H}^+]$$

The **Henderson-Hasselbalch** equation provides a convenient way to think about buffers and pH:

$$\text{pH} = \text{pKa} + \log \frac{[\text{A}^-]}{[\text{HA}]}$$

If one were examining the dissociation of acetic acid, the Henderson-Hasselbalch relationship, then, would be:

$$\text{pH} = \text{pKa} + \log \frac{[\text{CH}_3\text{COO}^-]}{[\text{CH}_3\text{COOH}]}$$

The Henderson-Hasselbalch equation can be used to determine if an aqueous solution of a conjugate acid/base pair is functioning as a buffer. If the concentration of the weak acid is equal to that of its conjugate base, the ratio of these two components is one. When this is the case, the Henderson-Hasselbalch equation reduces to

$$\text{pH} = \text{pKa}$$

because the  $\log(1)$  is equal to zero.

When the pH of the solution is equal to the pKa of the ionizing group, the solution is functioning at maximum buffering capacity (*best buffer*). An aqueous solution of a conjugate acid/base pair functions as a *good buffer* when the ratio of the conjugate base to weak acid ranges from 1:9 to 9:1. Substituting these ratios into the Henderson-Hasselbalch equation, one finds that this aqueous solution functions as a good buffer when the pH of the solution is within approximately one pH unit of the ionizing group's pKa.

$$\text{pH} = \text{pKa} \pm 1$$

because the  $\log(1/9)$  is -0.999 and the  $\log(9/1)$  is +0.999.

### Using the Henderson-Hasselbalch Equation

**Example A:** *This is a two-component buffer system meaning that the weak acid and its conjugate base are added separately.* How would you prepare 10mL of a 0.01M phosphate buffer, pH 7.40, from stock solutions of 0.10M  $\text{KH}_2\text{PO}_4$  and 0.25M  $\text{K}_2\text{HPO}_4$ ?  $\text{pKa}$  of  $\text{KH}_2\text{PO}_4 = 7.20$ . This problem is similar to Problem 1 in **Experimental Procedures**. The following approach may be helpful in solving this type of buffer problem in which both the conjugate acid and base are added separately. Please note that the numbers are not rounded off until the very end and are rounded based on the limits of the pipets, cylinders, etc. that are required to accurately measure the calculated volumes.

1. Use the Henderson Hasselbalch equation to find the ratio of  $\text{A}^-$  to HA.

$$\text{pH} = \text{pKa} + \log [\text{A}^-] / [\text{HA}]$$

$$7.40 = 7.20 + \log [\text{A}^-] / [\text{HA}]$$

$$0.20 = \log [\text{A}^-] / [\text{HA}]$$

$$1.584893192 = [\text{A}^-] / [\text{HA}]^*$$

\*Since  $[\text{A}^-] / [\text{HA}] = 1.584893192$ , we can say that  $[\text{A}^-] / [\text{HA}] = 1.584893192 / 1$ . In this case  $[\text{A}^-] = 1.584893192$ ;  $[\text{HA}] = 1$ .

2. Calculate the decimal fraction (part/whole) of each buffer component.

$$\text{A}^- = 1.584893192 / (1.000 + 1.584893192) = 1.584893192 / 2.584893192 = 0.61313682$$

$$\text{HA} = 1.000 / 2.584893192 = 0.38686318$$

3. Find the **molarity (M)** of each component in the buffer by simply multiplying the molarity of the buffer by the decimal fraction of each component.

$$M_{\text{A}^-} = 0.01\text{M} \times 0.61313682 = 0.006131368\text{M}$$

$$M_{\text{HA}} = 0.01\text{M} \times 0.38686318 = 0.003868632\text{M}$$

4. Calculate the moles of each component in the buffer. **Moles = Molarity x Liters** of buffer

$$\text{moles}_{\text{A}^-} = 0.006131368\text{M} \times 0.01\text{L} = 6.131 \times 10^{-5} \text{ moles}$$

$$\text{moles}_{\text{HA}} = 0.003868632\text{M} \times 0.01\text{L} = 3.869 \times 10^{-5} \text{ moles}$$

- Calculate the volume of each stock solution required to make the buffer

Liters of stock = moles of the buffer component / Molarity of the stock

$$L_{\text{A}^-} = 6.131 \times 10^{-5} \text{ moles} / 0.25 \text{ M} = 2.452 \times 10^{-4} \text{ L} = 245\mu\text{L}$$

$$L_{\text{HA}} = 3.869 \times 10^{-5} \text{ moles} / 0.10 \text{ M} = 3.869 \times 10^{-4} \text{ L} = 387\mu\text{L}$$

- To prepare this buffer, one would use appropriately-sized pipets to measure and transfer each component to a 10mL volumetric flask and bring the solution to volume with dH<sub>2</sub>O.

Example B: *This is a single-component buffer system. This means that one component of the conjugate acid/base pair will be generated in situ (in solution). In the following example, the conjugate base will be generated by reaction between the weak acid (acetic acid) and strong base (sodium hydroxide). How would you prepare 10mL of a 0.02M acetate buffer, pH 4.30, from stock solutions of 0.05M acetic acid (HAc) and 0.05M NaOH?  $pK_{\text{a}}_{\text{acetic acid}} = 4.76$ . This problem is similar to Problem 2 in **Experimental Procedures**.*

- Use the Henderson Hasselbalch equation to find the ratio of A<sup>-</sup> to HA.

$$\text{pH} = \text{pKa} + \log [\text{A}^-] / [\text{HA}]$$

$$4.30 = 4.76 + \log [\text{A}^-] / [\text{HA}]$$

$$-0.46 = \log [\text{A}^-] / [\text{HA}]$$

$$0.34673685 = [\text{A}^-] / [\text{HA}]$$

- Calculate the decimal fraction (part/whole) of each buffer component.

$$\text{A}^- = 0.34673685 / (1.00 + 0.34673685) = 0.34673685 / 1.34673685 = 0.257464441$$

$$\text{HA} = 1.00 / 1.34673685 = 0.742535559$$

- Find the **molarity (M)** of each component in the buffer by simply multiplying the molarity of the buffer by the decimal fraction of each component.

$$M_{\text{A}^-} = 0.02\text{M} \times 0.257464441 = 0.005149289 = 5.15 \times 10^{-3} \text{ M}$$

$$M_{\text{HA}} = 0.02\text{M} \times 0.742535559 = 0.014850711 = 1.49 \times 10^{-2} \text{ M}$$

- Calculate the moles of each component in the buffer. Moles = Molarity x Liters of buffer

$$\text{moles}_{\text{A}^-} = 0.005149289\text{M} \times 0.01\text{L} = 5.14929 \times 10^{-5} \text{ moles}$$

$$\text{moles}_{\text{HA}} = 0.014850711\text{M} \times 0.01\text{L} = 1.48507147 \times 10^{-4} \text{ moles}$$

- Note: Since this buffer is prepared by the reaction of a weak acid (HAc) with a strong base (NaOH), you must determine the **total moles** of the weak acid component needed because the conjugate base is made *in situ*.

$$\text{Total moles} = 5.15 \times 10^{-5} \text{ moles NaOH} + 1.49 \times 10^{-4} \text{ moles HAc} = \mathbf{2.00 \times 10^{-4}}$$

**moles HAc.** *This sum indicates that, although in the buffer one only needs  $1.49 \times 10^{-4}$  moles HA, an additional  $5.15 \times 10^{-5}$  moles is needed to generate the conjugate base in situ.*

- Calculate the volume of each stock solution required to make the buffer

Liters of stock = Moles of the buffer component / Molarity of the stock

$$L_{A^-} = 5.14929 \times 10^{-5} \text{ moles} / 0.05 \text{ M} = 1.03 \times 10^{-3} \text{ L NaOH} = 1.0 \text{ mL}$$

$$L_{\text{HA}} = \mathbf{2.00 \times 10^{-4} \text{ moles}} / 0.05 \text{ M} = 4.00 \times 10^{-3} \text{ L CH}_3\text{COOH} = 4.0 \text{ mL}$$

7. To prepare this buffer, one would use appropriately sized pipets or cylinders to measure and transfer each component to a 10 mL volumetric flask and bring the solution to volume with dH<sub>2</sub>O.


**Example C:** What is the pH of a solution resulting from a mixture of 200 mL of 0.1 M NaOH and 100 mL of 0.3 M HAc? This problem is similar to Problem 3 in **Experimental Procedures**. It is a *limiting reagent* type problem.

1. Calculate the initial moles of each reactant. Moles = Molarity x Liters

$$\text{Moles}_{\text{NaOH}} = 0.1 \text{ M} \times 0.2 \text{ L} = 0.02 \text{ moles}$$

$$\text{Moles}_{\text{HAc}} = 0.3 \text{ M} \times 0.1 \text{ L} = 0.03 \text{ moles}$$

2. Write the balanced equation for this reaction. *Note: This is a neutralization reaction.*



3. Determine the moles of each reactant remaining after the reaction. *Recall: This is a limiting reagent type problem.* Since NaOH is present in the least amount and, according to the balanced equation the reactants react in a 1:1 ratio, the NaOH is completely consumed (i.e., zero moles of NaOH remain). Because the reactants react in a 1:1 ratio,

$$\begin{aligned} \text{Moles}_{\text{HAc remaining}} &= \text{Moles}_{\text{HAc initial}} - \text{Moles}_{\text{HAc reacted}} \\ &= 0.03 \text{ moles} - 0.02 \text{ moles} = 0.01 \text{ moles} \end{aligned}$$

4. Determine the molarity of the buffer components (i.e., the conjugate acid-base pair).

$$M_{\text{HAc}} = \text{moles}_{\text{HAc remaining}} / \text{Liters solution} = 0.01 \text{ moles} / 0.3 \text{ L} = 0.033 \text{ M}$$

The molarity of the conjugate base Ac<sup>-</sup> resulting from the dissolution of the product NaAc is:

$$M_{\text{Ac}^-} = \text{moles formed} / \text{Liters solution} = 0.02 \text{ moles} / 0.3 \text{ L} = 0.067 \text{ M}$$

5. Use the H-H equation to calculate the pH.

$$\begin{aligned} \text{pH} &= \text{pK}_a + \log [A^-] / [HA] \\ &= 4.76 + \log (0.067/0.033) \\ &= 4.76 + \log 2.02 \\ &= 4.76 + 0.31 \\ &= 5.07 \end{aligned}$$

**Try this one: Determine the pH of a solution resulting from a mixture of 5mL of 0.1M KOH and 30mL of 0.1M HAc.**

### Experimental Procedures:

1. Prepare 10 mL of a 0.01 M phosphate buffer, pH 7.70, from stock solutions of 0.1 M  $\text{K}_2\text{HPO}_4$  and 0.2 M  $\text{KH}_2\text{PO}_4$ . ( $\text{pK}_a$  for the weak acid = 7.20).

- A. Use the Henderson-Hasselbalch equation to calculate the volume of each stock solution needed.

$$\text{pH} = \text{pK}_a + \log \left[ \frac{\text{conjugate base}}{\text{weak acid}} \right]$$

- B. Check your calculations with other students. See the instructor if there is uncertainty.
- C. Make the solution and check the pH of a portion of your buffer solution using the pH meter. (Instructions are on **page 33**)

2. Prepare 10 mL of 0.01 M acetate buffer, pH 3.80, from stock solutions of 0.1 M acetic acid and 0.02 M sodium hydroxide.  $\text{pK}_{a_{\text{acetic acid}}} = 4.76$ .

- A. – C. See above.

- D. Calculate the exact volume of the 0.01 M acetate buffer required to make 10 mL of a 0.0005 M acetate buffer.

1. Prepare this new buffer using the following equation to aid you in your calculations.

$$V_1 = \frac{(V_2 \times M_2)}{M_1}$$

where:  $V_1$  = the volume of the concentrated solution (liters, L)

$V_2$  = the volume of the diluted solution (liters, L)


$M_1$  = the molar concentration of the concentrated solution (moles/L)

$M_2$  = the molar concentration of the diluted solution (moles/L)

2. Check the pH of this new buffer.

- E. Calculate the exact volume of the 0.01 M acetate buffer required to make 10 mL of a 0.001 M acetate buffer.

1. Prepare this new buffer.
2. Check its pH.



3. What is the pH of a solution resulting from the addition of 2.5 mL of 0.02M NaOH to 5.0 mL 0.25M acetic acid? Use the Henderson-Hasselbalch equation to calculate the pH of the resulting solution, and then follow directions B and C above.  $pK_{\text{acetic acid}} = 4.76$ .

**Record the data in your notebook (you should have 5 pH values):**

## The Preparation of Buffer Solutions of Desired pH STUDY GUIDE

1. The hydrogen ion concentration of tomato juice (pH 4) is how many times higher than the hydrogen ion concentration of 0.01 M NaOH (pH 12)?
2. Phosphoric acid  $\text{H}_3\text{PO}_4$  has three dissociable protons (it is triprotic) thus three pKa values; 2.12, 7.21, and 12.32). What volume of 200 mM phosphoric acid and what volume of 500 mM NaOH are required to make 100 mL of 20 mM phosphate buffer pH 12?
3. What must be done to the electrode of the pH meter before placing it into a solution to measure its pH?
4. In the phosphate buffer system containing  $\text{K}_2\text{HPO}_4$  and  $\text{KH}_2\text{PO}_4$ , what is the weak acid? What is its conjugate base?
5. If the weak acid in a buffer system has pKa 10.5, in what pH range is the system a most effective buffer?
6. On the laboratory shelf are 250mM solutions of both acetic acid and sodium hydroxide. How would you make a 100 mL solution of 25mM acetate buffer of pH 5.50 using these stock solutions?
7. How would you prepare a liter of 25% NaCl (w:w) aqueous solution?
8. How would you prepare a 1.5% NaBr (w:v) aqueous solution?
9. How would you prepare a 65% acetone (v:v) aqueous solution?
10. How would you prepare 500 mL of 0.12 $\mu\text{M}$  aqueous solution of glucose ( $\text{C}_6\text{H}_{12}\text{O}_6$ )?
11. What is the percent conjugate acid in an acetate buffer of pH 3.80? (pKa = 4.76)
12. How many  $\mu\text{moles}$  are equivalent to 100 pmoles?
- 13.
14. What is the ratio of conjugate base to weak acid in a buffer when the pH drops to two units below the pKa?