HEWLETT-PACKARD

NEWS FOR
USERS OF
ANALYTICAL
INSTRUMENTATION

IN THIS ISSUE

- ☐ Why Skunks Smell So Bad
- Screening Samples With SFC
 - New Ways for Automated On-Column Injection
- ☐ Customizing LC for On-Line Quality Control
- ☐ Searching for a Dual-Purpose Column
- ☐ 3-D HPLC for Complex Samples

IDENTIFICATION OF THE VOLATILE COMPONENTS IN SKUNK SPRAY BY GC

Why Skunks Smell So Bad

by William F. Wood, Department of Chemistry, Humboldt State University, Arcata (California)

Natural product chemists have long been interested in the constituents of the defensive spray of the striped skunk, *Mephitis mephitis*. GC/MS analysis of this repulsive-smelling secretion shows it to contain seven major components, including thiols, alkaloids, and thioacetate derivatives of the thiols.


History of Skunk Research

Research on the repulsive-smelling defensive secretion of the striped skunk has long been of interest to chemists. In 1862, the German chemist Swarts reported the secretion to be composed of 18% sulfur, to have several fractions that distilled at different temperatures and contain a nitrous base. In 1897, 35 years later, Aldrich and Jones of Johns Hopkins University showed this base to be 2-methylquinoline (Compound 5 — see table). 2

A previous publication by Aldrich in 1896 incorrectly identified 1-butanethiol (butyl mercaptan) as the major constituent in this secretion.³ This error was widely perpetuated for over 80 years, until 1975, when Andersen and Bernstein showed the compound to be *trans*-2-butene-1-thiol (Compound 1).⁴ In the latter investigation, a second major odoriferous compound, 3-methyl-1-butanethiol (Compound 2) was identified.

The present study of the defensive secretion of the striped skunk has shown that there are seven major components in this secretion (see figure). Besides the thiols identified by Andersen and Bernstein and the 2-methylquinoline found by Aldrich, the secretion was shown to contain a number of other compounds. These include the alkaloid 2-quinolinemethanethiol and thioacetate derivatives of all thiols.

Major components of the striped skunk's defensive spray in order of retention time.

Capillary GC chromatogram of striped-skunk defensive secretion (numerals represent Compound No.).

Collection of Secretion

The skunks used in this research were live-trapped in areas where they were a nuisance. After treatment with ether, a blunt needle was gently inserted into the anal sac, and the contents were removed by suction. A few drops of the secretion were diluted with 50 ml of diethyl ether, and 2 μ l of this solution were immediately analyzed by GC/MS.

Analysis of the Sample

The capillary GC chromatogram in the figure was obtained by analysis of the ether solution using an HP 5890 gas chromatograph fitted with an HP 5970 MSD. Integration of the TIC showed seven compounds to be present at a concentration greater than 1%. The variation of the amounts of each component from four different skunks is shown in the table.

Identification of Components

The first component, trans-2-butene-1-thiol, was identified by its mass spectrum. The second compound, 3-methyl-1-butanethiol, was identified by comparison with an entry in the

NBS Mass Spectral Library (Revision E). It was confirmed by comparison with an authentic sample of this compound. The third compound, trans-2-butenyl thioacetate (Compound 3) was identified by comparison with a synthetically prepared sample. Compound 4, 3-methylbutenyl thioacetate, was likewise identified.

2-Methylquinoline had mass spectra and retention time identical to those of an authentic sample. The two other alkaloids, 2-quinolinemethanethiol (Compound 6) and 2-quinolinemethyl thioacetate (Compound 7) were identified by comparison with synthetically prepared samples.

Why Skunk Odor Persists

"Skunk odor" can linger for days on sprayed animals and will return weeks later when the animal's fur is wet. The persistence of the odor may be attributed to the thioacetates in the secretion. These compounds decompose slowly in water and continue

The Instrumentation Used

HP 5890 gas chromatograph Column: HP-1, 12 m x 0.2 mm x 0.33 μm Temperature: 400°C for 4 min, 40° to 250°C at 25°C/min, 250°C 4 min

HP 5970 mass spectral detector HP 9133 ChemStation with NBS Mass Spectral Library (Revision E)

Percent Variation of Components (4 Skunks)		
Compound No.	Amount (% TIC)	
1	38-40	
2	18-26	
3	12-18	
4	2-3	
5	4-11	
6	4-12	
7	14	

to release the more potent thiols for days. To speed the destruction of the thioacetates, animals can be washed with laundry or castile soap. Thioacetates are rapidly hydrolyzed under mildly basic conditions.

About the Author

William Wood is a professor of chemistry at Humboldt State University. His research interests are centered on the chemistry of biologiin the field of semiochemistry. His published works in this field include the identification of arthropod secretions), the investigation of mammalian defensive secretions, and the identification of natural products from plants. He was a research fellow at the International Centre for Insect Physiology and Ecology (ICIPE) in Nairobi, Kenya.

References

- 1. Swarts, Justus Liebigs Ann. Chem. 123, 266
- 2. Aldrich, T. B., and Jones, W., J. Exp. Med. 2, 439 (1897).
- 3. Aldrich, T. B., J. Exp. Med. 1, 323 (1896).
- 4. Andersen, K. K. and Bernstein, D. T., J. Chem. Ecol. 1, 493 (1975).
- 5. Wood, W. F., J. Chem. Ecol. 16, 2057 (1990)

cally important compounds, primarily pheromones and allomones (defensive

More Power for HPLC ChemStation

Now that the new HP 9000 Series 300 Model 362 computer is being used in the HPLC ChemStation (Pascal series), performance of the HP 1090 Series II/M HPLC system has received still another boost. The computer, running at 25 MHz and offering direct access to the built-in hard drive, allows data evaluation virtually simultaneously with data acquisition.

The ChemStation will be offered with a VGA color monitor, 4 Mbytes of computer RAM, a built-in 3.5-inch 1.44-Mbyte flexible-disk drive, and a builtin 200-Mbyte hard-disk drive.

For more information on the new HPLC ChemStation, check A on the PEAK reply card.

Why Skunks Smell So Bad **Screening Samples With SFC** 5 **New Savings With UV-Vis** Solutions **New Ways for Automated On-Column Injection** The Best Benchtop GC/MS Just 7 **Got Better**

10 Searching for a Dual-Purpose Column

Customizing LC for On-Line

Quality Control

8

11 3-D HPLC for Complex Samples

12 **New Analytical Tools from HP**

PEAK is published quarterly by the Hewlett-Packard Company and is distributed to users of analytical instrumentation throughout the world.

If you are not already receiving your own copy of PEAK, please check the appropriate box on the reply card.

Editor: Hans Toepfer Hewlett-Packard Company **Analytical Products Group** 3495 Deer Creek Road Palo Alto, CA 94304 U.S.A.

Cover Illustration: Dick Cole Graphic Design: Guild West, Los Altos, California, U.S.A. Copyright o Hewlett-Packard Company 1992

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